



## User's Manual

# Human sHLA-G ELISA Kit



DEIA2208



96T



This product is for research use only and is not intended for diagnostic use.

For illustrative purposes only. To perform the assay the instructions for use provided with the kit have to be used.

---

### Creative Diagnostics

 Address: 45-1 Ramsey Road, Shirley, NY 11967, USA

 Tel: 1-631-624-4882 (USA) 44-161-818-6441 (Europe)  Fax: 1-631-938-8221

 Email: [info@creative-diagnostics.com](mailto:info@creative-diagnostics.com)  Web: [www.creative-diagnostics.com](http://www.creative-diagnostics.com)

---

## PRODUCT INFORMATION

### Intended Use

The sHLA-G ELISA is a sandwich enzyme immunoassay for the measurement of soluble forms of human leukocyte antigen-G ( sHLA-G ).

### Storage

Store all contents at 2-8°C. For more detailed information, please download the following document on our website.

### Detection Range

3.91-125 U/mL

### Sensitivity

3 U/mL

