



User's Manual

GLP-1 Total (Multi-Species) ELISA Kit



DEIA2216



96T






This product is for research use only and is not intended for diagnostic use.

For illustrative purposes only. To perform the assay the instructions for use provided with the kit have to be used.

Creative Diagnostics

 Address: 45-1 Ramsey Road, Shirley, NY 11967, USA

 Tel: 1-631-624-4882 (USA) 44-161-818-6441 (Europe)  Fax: 1-631-938-8221

 Email: info@creative-diagnostics.com  Web: www.creative-diagnostics.com

PRODUCT INFORMATION

Intended Use

This GLP-1 Total ELISA kit is used for the non-radioactive quantification of GLP-1 (7-36 and 9-36) in serum, plasma, and cell tissue culture.

Storage

Store all contents at 2-8°C. For more detailed information, please download the following document on our website.

Detection Range

4.1-1000 pM

Sensitivity

1.5 pmol/L

