



User's Manual

Trypsin ELISA Kit



DEIA2227



192T



This product is for research use only and is not intended for diagnostic use.

For illustrative purposes only. To perform the assay the instructions for use provided with the kit have to be used.

Creative Diagnostics

 Address: 45-1 Ramsey Road, Shirley, NY 11967, USA

 Tel: 1-631-624-4882 (USA) 44-161-818-6441 (Europe)  Fax: 1-631-938-8221

 Email: info@creative-diagnostics.com  Web: www.creative-diagnostics.com

PRODUCT INFORMATION

Intended Use

The Neonate Trypsin ELISA kit is for measurement of human immunoreactive trypsin (IRT) from blood spot samples collected on Schleicher and Schuell's filter paper.

Storage

When stored at 2-8°C unopened reagents will retain reactivity until expiration date. For more detailed information, please download the following document on our website.