



**User's Manual**

# Human FSH (Follicle Stimulating Hormone) ELISA Kit

REF

DEIA2229



96T

RUO

This product is for research use only and is not intended for diagnostic use.

For illustrative purposes only. To perform the assay the instructions for use provided with the kit have to be used.

---

**Creative Diagnostics**

 **Address:** 45-1 Ramsey Road, Shirley, NY 11967, USA

 **Tel:** 1-631-624-4882 (USA) 44-161-818-6441 (Europe)  **Fax:** 1-631-938-8221

 **Email:** [info@creative-diagnostics.com](mailto:info@creative-diagnostics.com)  **Web:** [www.creative-diagnostics.com](http://www.creative-diagnostics.com)

---

## PRODUCT INFORMATION

### Intended Use

The FSH ELISA is an enzyme immunoassay for the quantitative in vitro diagnostic measurement of Follicle Stimulating Hormone (FSH) in serum.

### Storage

Store the kit at 2-8°C. Keep microwells sealed in a dry bag with desiccants. For more detailed information, please download the following document on our website.

### Detection Range

0.86-100 mIU/mL

### Sensitivity

0.856 mIU/mL