



User's Manual

Serotonin ELISA Kit



DEIA2236



96T



This product is for research use only and is not intended for diagnostic use.

For illustrative purposes only. To perform the assay the instructions for use provided with the kit have to be used.

Creative Diagnostics

 Address: 45-1 Ramsey Road, Shirley, NY 11967, USA

 Tel: 1-631-624-4882 (USA) 44-161-818-6441 (Europe)  Fax: 1-631-938-8221

 Email: info@creative-diagnostics.com  Web: www.creative-diagnostics.com

PRODUCT INFORMATION

Intended Use

Enzyme immunoassay for the in-vitro diagnostic quantitative determination of Serotonin in human serum, plasma, platelets and urine. Further the test can be used for research of tissue homogenates and cell culture supernatants.

Storage

The kit is shipped at ambient temperature and should be stored at 2-8°C. For more detailed information, please download the following document on our website.

Sensitivity

0.014 ng/mL



Tel: 1-631-624-4882 (USA)

Tel: 44-161-818-6441 (Europe)



Fax: 1-631-938-8221



Email: info@creative-diagnostics.com