



User's Manual

Human alpha-Atrial Natriuretic Polypeptide ELISA Kit

REF

DEIA2244



96T

RUO

This product is for research use only and is not intended for diagnostic use.

For illustrative purposes only. To perform the assay the instructions for use provided with the kit have to be used.

Creative Diagnostics

 **Address: 45-1 Ramsey Road, Shirley, NY 11967, USA**

 **Tel: 1-631-624-4882 (USA) 44-161-818-6441 (Europe)**  **Fax: 1-631-938-8221**

 **Email: info@creative-diagnostics.com**  **Web: www.creative-diagnostics.com**

PRODUCT INFORMATION

Intended Use

This Enzyme Immunoassay kit is designed to detect a specific peptide and its related peptides based on the principle of competitive enzyme immunoassay.

Storage

Store all contents at 4°C. For more detailed information, please download the following document on our website.