



User's Manual

Food Allergies IgG - Allerquant 90G ELISA Kit



DEIA2264



3 x 96T



This product is for research use only and is not intended for diagnostic use.

For illustrative purposes only. To perform the assay the instructions for use provided with the kit have to be used.

Creative Diagnostics

 **Address: 45-1 Ramsey Road, Shirley, NY 11967, USA**

 **Tel: 1-631-624-4882 (USA) 44-161-818-6441 (Europe)**  **Fax: 1-631-938-8221**

 **Email: info@creative-diagnostics.com**  **Web: www.creative-diagnostics.com**

PRODUCT INFORMATION

Intended Use

The 90 Food IgG Elisa Test is for measuring the relative amount of food-specific IgG antibody in human serum.

Storage

For more detailed information, please download the following document on our website.

Detection Range

50-400 U/mL