



## User's Manual

# Human CMV IgM ELISA Kit



DEIA2274



96T



This product is for research use only and is not intended for diagnostic use.

For illustrative purposes only. To perform the assay the instructions for use provided with the kit have to be used.

---

### Creative Diagnostics

 Address: 45-1 Ramsey Road, Shirley, NY 11967, USA

 Tel: 1-631-624-4882 (USA) 44-161-818-6441 (Europe)  Fax: 1-631-938-8221

 Email: [info@creative-diagnostics.com](mailto:info@creative-diagnostics.com)  Web: [www.creative-diagnostics.com](http://www.creative-diagnostics.com)

---

## PRODUCT INFORMATION

### Intended Use

Enzyme Immunoassay for the Detection of IgM antibodies to Cytomegalovirus in Human Serum.

### Storage

Store the kit at 2-8°C. Keep microwells sealed in a dry bag with desiccants. For more detailed information, please download the following document on our website.