



User's Manual

Human 'TGF beta 2 ELISA Kit

REF DEIA2303

Σ 96T

RUO

This product is for research use only and is not intended for diagnostic use.

For illustrative purposes only. To perform the assay the instructions for use provided with the kit have to be used.

Creative Diagnostics

 Address: 45-1 Ramsey Road, Shirley, NY 11967, USA

 Tel: 1-631-624-4882 (USA) 44-161-818-6441 (Europe)  Fax: 1-631-938-8221

 Email: info@creative-diagnostics.com  Web: www.creative-diagnostics.com

PRODUCT INFORMATION

Intended Use

The TGF- β 2 ELISA is an enzyme immunoassay for the quantitative in vitro diagnostic measurement of TGF- β 2 in serum, plasma and cell culture supernatant.

Storage

Store the reagents at 2-8°C until expiration date. For more detailed information, please download the following document on our website.

Detection Range

10-1000 pg/mL

Sensitivity

10 pg/mL