



User's Manual

Human ANA IgG Screen ELISA



DEIA2356



96T



This product is for research use only and is not intended for diagnostic use.

For illustrative purposes only. To perform the assay the instructions for use provided with the kit have to be used.

Creative Diagnostics

 Address: 45-1 Ramsey Road, Shirley, NY 11967, USA

 Tel: 1-631-624-4882 (USA) 44-161-818-6441 (Europe)  Fax: 1-631-938-8221

 Email: info@creative-diagnostics.com  Web: www.creative-diagnostics.com

PRODUCT INFORMATION

Intended Use

This ELISA test system is designed to detect IgG class antibodies to a variety of common nuclear antigens in human sera.

Storage

The kit should be stored at 2-8°C. For more detailed information, please download the following document on our website.

Sensitivity

Sensitivity is 90.4%.