



User's Manual

Cardiolipin IgG/IgA/IgM ELISA Kit



DEIA2360



96T



This product is for research use only and is not intended for diagnostic use.

For illustrative purposes only. To perform the assay the instructions for use provided with the kit have to be used.

Creative Diagnostics

 **Address: 45-1 Ramsey Road, Shirley, NY 11967, USA**

 **Tel: 1-631-624-4882 (USA) 44-161-818-6441 (Europe)**  **Fax: 1-631-938-8221**

 **Email: info@creative-diagnostics.com**  **Web: www.creative-diagnostics.com**

PRODUCT INFORMATION

Intended Use

This product is intended for the detection and qualitative determination of IgG, IgA, and IgM antibodies to cardiolipin in human sera.

Storage

Store unopened kit between 2 and 8°C. For more detailed information, please download the following document on our website.

Sensitivity

Relative Sensitivity = 35/35 = 100%