



User's Manual

Human Allergen Specific IgE ELISA Kit



DEIA2429



2 x 96T



This product is for research use only and is not intended for diagnostic use.

For illustrative purposes only. To perform the assay the instructions for use provided with the kit have to be used.

Creative Diagnostics

 **Address: 45-1 Ramsey Road, Shirley, NY 11967, USA**

 **Tel: 1-631-624-4882 (USA) 44-161-818-6441 (Europe)**  **Fax: 1-631-938-8221**

 **Email: info@creative-diagnostics.com**  **Web: www.creative-diagnostics.com**

PRODUCT INFORMATION

Intended Use

To qualitate and quantitate allergen specific human Immunoglobulin E.

Storage

Store all components at 2-8°C. Do not freeze all or in part. For more detailed information, please download the following document on our website.

Detection Range

0.35-0.70 PRU/mL, depending on the type of allergen used.