



User's Manual

Human Neo-Natal TSH ELISA Kit



DEIA2475



96T



This product is for research use only and is not intended for diagnostic use.

For illustrative purposes only. To perform the assay the instructions for use provided with the kit have to be used.

Creative Diagnostics

 **Address: 45-1 Ramsey Road, Shirley, NY 11967, USA**

 **Tel: 1-631-624-4882 (USA) 44-161-818-6441 (Europe)**  **Fax: 1-631-938-8221**

 **Email: info@creative-diagnostics.com**  **Web: www.creative-diagnostics.com**

PRODUCT INFORMATION

Intended Use

The Quantitative Determination of Thyrotropin Concentration in Human Whole Blood by a Microplate Immunoassay.

Storage

All reagents provided are stored at 4°C. For more detailed information, please download the following document on our website.

Detection Range

Range 0.7 µIU/mL-34 µIU/mL

Sensitivity

1.0 µIU/mL



Tel: 1-631-624-4882 (USA)

Tel: 44-161-818-6441 (Europe)



Fax: 1-631-938-8221



Email: info@creative-diagnostics.com