



User's Manual

Campylobacter ELISA Kit (Fecal)



DEIA2560



96T



This product is for research use only and is not intended for diagnostic use.

For illustrative purposes only. To perform the assay the instructions for use provided with the kit have to be used.

Creative Diagnostics

 **Address: 45-1 Ramsey Road, Shirley, NY 11967, USA**

 **Tel: 1-631-624-4882 (USA) 44-161-818-6441 (Europe)**  **Fax: 1-631-938-8221**

 **Email: info@creative-diagnostics.com**  **Web: www.creative-diagnostics.com**

PRODUCT INFORMATION

Intended Use

This microwell enzyme linked immunoabsorbant assay (ELISA) detection kit is a diagnostic immunoassay for the detection of antigen to Campylobacter species in human feces using peroxidase as the indicator enzyme.

Storage

Reagents, strips and bottled components: Store at 2-8°C. Squeeze bottle containing diluted wash buffer may be stored at room temperature. For more detailed information, please download the following document on our website.

Sensitivity

This assay can detect approximately 104 to 105 CFU per mL of feces.