



User's Manual

E.coli Antigen In Food ELISA Kit



DEIA2562



96T



This product is for research use only and is not intended for diagnostic use.

For illustrative purposes only. To perform the assay the instructions for use provided with the kit have to be used.

Creative Diagnostics

 **Address: 45-1 Ramsey Road, Shirley, NY 11967, USA**

 **Tel: 1-631-624-4882 (USA) 44-161-818-6441 (Europe)**  **Fax: 1-631-938-8221**

 **Email: info@creative-diagnostics.com**  **Web: www.creative-diagnostics.com**

PRODUCT INFORMATION

Intended Use

This assay is an enzyme-linked immunosorbent assay (ELISA) that may be used to screen food products for the presence of E. coli O157 antigen.

Storage

Reagents, strips and bottled components: Store at 2-8°C. Squeeze bottle containing diluted wash buffer may be stored at room temperature. For more detailed information, please download the following document on our website.