



User's Manual

Mouse /Rat/Porcine/Canine TGF- β 1 ELISA Kit



DEIA2657



96T



This product is for research use only and is not intended for diagnostic use.

For illustrative purposes only. To perform the assay the instructions for use provided with the kit have to be used.

Creative Diagnostics

 **Address: 45-1 Ramsey Road, Shirley, NY 11967, USA**

 **Tel: 1-631-624-4882 (USA) 44-161-818-6441 (Europe)**  **Fax: 1-631-938-8221**

 **Email: info@creative-diagnostics.com**  **Web: www.creative-diagnostics.com**

PRODUCT INFORMATION

Intended Use

For the quantitative determination of activated Mouse /rat/porcine/canine Transforming Growth Factor beta 1 (TGF- β 1) concentrations in cell culture supernates, serum, and plasma.

Storage

Store the reagents at 2-8°C until expiration date. For more detailed information, please download the following document on our website.