



**User's Manual**

# Human sTNF RII/TNFRSF1B ELISA Kit



**DEIA2695**



**96T**



This product is for research use only and is not intended for diagnostic use.

For illustrative purposes only. To perform the assay the instructions for use provided with the kit have to be used.

---

**Creative Diagnostics**

 **Address: 45-1 Ramsey Road, Shirley, NY 11967, USA**

 **Tel: 1-631-624-4882 (USA) 44-161-818-6441 (Europe)**  **Fax: 1-631-938-8221**

 **Email: [info@creative-diagnostics.com](mailto:info@creative-diagnostics.com)**  **Web: [www.creative-diagnostics.com](http://www.creative-diagnostics.com)**

---

## PRODUCT INFORMATION

### Intended Use

For the quantitative determination of Human soluble tumor necrosis factor receptor II (sTNF RII) concentrations in cell culture supernates, serum, plasma, and urine.

### Storage

Store the reagents at 2-8°C until expiration date. For more detailed information, please download the following document on our website.

### Detection Range

0.2-2.3 pg/mL

### Sensitivity

0.6 pg/mL



**Tel:** 1-631-624-4882 (USA)

**Tel:** 44-161-818-6441 (Europe)



**Fax:** 1-631-938-8221



**Email:** [info@creative-diagnostics.com](mailto:info@creative-diagnostics.com)