



## User's Manual

# Human PDGF-BB ELISA Kit



DEIA2714



96T



This product is for research use only and is not intended for diagnostic use.

For illustrative purposes only. To perform the assay the instructions for use provided with the kit have to be used.

---

### Creative Diagnostics

 Address: 45-1 Ramsey Road, Shirley, NY 11967, USA

 Tel: 1-631-624-4882 (USA) 44-161-818-6441 (Europe)  Fax: 1-631-938-8221

 Email: [info@creative-diagnostics.com](mailto:info@creative-diagnostics.com)  Web: [www.creative-diagnostics.com](http://www.creative-diagnostics.com)

---

## PRODUCT INFORMATION

### Intended Use

For the quantitative determination of human Platelet-Derived Growth Factor-BB (PDGF-BB) concentrations in cell culture supernates, serum, and platelet-poor plasma.

### Storage

Store at 2-8°C. Do not use past kit expiration date. For more detailed information, please download the following document on our website.

### Sensitivity

15 pg/mL



**Tel:** 1-631-624-4882 (USA)

**Tel:** 44-161-818-6441 (Europe)



**Fax:** 1-631-938-8221



**Email:** [info@creative-diagnostics.com](mailto:info@creative-diagnostics.com)