



User's Manual

Human Soluble VEGF R2 ELISA Kit



DEIA2723



96T



This product is for research use only and is not intended for diagnostic use.

For illustrative purposes only. To perform the assay the instructions for use provided with the kit have to be used.

Creative Diagnostics

 **Address: 45-1 Ramsey Road, Shirley, NY 11967, USA**

 **Tel: 1-631-624-4882 (USA) 44-161-818-6441 (Europe)**  **Fax: 1-631-938-8221**

 **Email: info@creative-diagnostics.com**  **Web: www.creative-diagnostics.com**

PRODUCT INFORMATION

Intended Use

For the quantitative determination of human soluble VEGF Receptor 2 (sVEGF R2) concentrations in cell culture supernates, cell lysates, serum, and plasma.

Storage

Store the unopened kit at 4°C upon receipt and when it is not in use. For more detailed information, please download the following document on our website.

Detection Range

1.0-11.4 pg/mL

Detection Limit

4.6 pg/mL

