



User's Manual

Mouse Lipocalin-2/NGAL ELISA Kit



DEIA2730



96T



This product is for research use only and is not intended for diagnostic use.

For illustrative purposes only. To perform the assay the instructions for use provided with the kit have to be used.

Creative Diagnostics

 Address: 45-1 Ramsey Road, Shirley, NY 11967, USA

 Tel: 1-631-624-4882 (USA) 44-161-818-6441 (Europe)  Fax: 1-631-938-8221

 Email: info@creative-diagnostics.com  Web: www.creative-diagnostics.com

PRODUCT INFORMATION

Intended Use

For the quantitative determination of mouse Lipocalin-2/NGAL concentrations in cell culture supernates, mouse serum, platelet-poor plasma, and urine.

Storage

Store at 2-8°C. Do not use past kit expiration date.

Precision

Intra-assay Precision: 4.1%-6.7%

Inter-assay Precision: 6.0%-7.2%

Detection Range

1.40-8.80 pg/mL

Sensitivity

3.18 pg/mL