



**User's Manual**

# Human TIM-1/KIM-1/HAVCR ELISA Kit



DEIA2760



96T



This product is for research use only and is not intended for diagnostic use.

For illustrative purposes only. To perform the assay the instructions for use provided with the kit have to be used.

---

**Creative Diagnostics**

 **Address: 45-1 Ramsey Road, Shirley, NY 11967, USA**

 **Tel: 1-631-624-4882 (USA) 44-161-818-6441 (Europe)**  **Fax: 1-631-938-8221**

 **Email: [info@creative-diagnostics.com](mailto:info@creative-diagnostics.com)**  **Web: [www.creative-diagnostics.com](http://www.creative-diagnostics.com)**

---

## PRODUCT INFORMATION

### Intended Use

For the quantitative determination of human T cell Immunoglobulin-Mucin (TIM-1) concentrations in urine.

### Storage

Store at 2-8°C. Do not use past kit expiration date. For more detailed information, please download the following document on our website.

### Detection Range

0.003-0.046 ng/mL

### Detection Limit

0.009 ng/mL.