



**User's Manual**

# Human CXCL13/BLC/BCA-1 ELISA Kit



DEIA2768



96T



This product is for research use only and is not intended for diagnostic use.

For illustrative purposes only. To perform the assay the instructions for use provided with the kit have to be used.

---

**Creative Diagnostics**

 **Address: 45-1 Ramsey Road, Shirley, NY 11967, USA**

 **Tel: 1-631-624-4882 (USA) 44-161-818-6441 (Europe)**  **Fax: 1-631-938-8221**

 **Email: [info@creative-diagnostics.com](mailto:info@creative-diagnostics.com)**  **Web: [www.creative-diagnostics.com](http://www.creative-diagnostics.com)**

---

## PRODUCT INFORMATION

### Intended Use

For the quantitative determination of Human B-lymphocyte Chemoattractant/B Cell-attracting Chemokine 1 (BLC/BCA-1) concentrations in cell culture supernates, serum, plasma, and saliva.

### Storage

Store the reagents at 2-8°C until expiration date. For more detailed information, please download the following document on our website.

### Detection Range

0.43-3.97 pg/mL

### Sensitivity

1.64 pg/mL