



User's Manual

# Human COMP ELISA Kit



DEIA2816



96T



This product is for research use only and is not intended for diagnostic use.

For illustrative purposes only. To perform the assay the instructions for use provided with the kit have to be used.

---

Creative Diagnostics

 Address: 45-1 Ramsey Road, Shirley, NY 11967, USA

 Tel: 1-631-624-4882 (USA) 44-161-818-6441 (Europe)  Fax: 1-631-938-8221

 Email: [info@creative-diagnostics.com](mailto:info@creative-diagnostics.com)  Web: [www.creative-diagnostics.com](http://www.creative-diagnostics.com)

---

## PRODUCT INFORMATION

### Intended Use

For the quantitative determination of human Cartilage Oligomeric Matrix Protein (COMP) concentrations in cell culture supernates, cell lysates, serum, and plasma.

### Storage

Store at 2-8°C. Do not use past kit expiration date. For more detailed information, please download the following document on our website.

### Detection Range

0.005-0.036 ng/mL

### Detection Limit

0.010 ng/mL