



User's Manual

Mouse sTNF RII/TNFRSF1B ELISA Kit



DEIA2845



96T



This product is for research use only and is not intended for diagnostic use.

For illustrative purposes only. To perform the assay the instructions for use provided with the kit have to be used.

Creative Diagnostics

 **Address: 45-1 Ramsey Road, Shirley, NY 11967, USA**

 **Tel: 1-631-624-4882 (USA) 44-161-818-6441 (Europe)**  **Fax: 1-631-938-8221**

 **Email: info@creative-diagnostics.com**  **Web: www.creative-diagnostics.com**

PRODUCT INFORMATION

Intended Use

For the quantitative determination of Mouse * soluble tumor necrosis factor receptor II (sTNF RII) concentrations in cell culture supernates, Mouse serum/plasma, and rat serum/plasma.

*This assay also recognizes rat sTNF RII.

Storage

Store the reagents at 2-8°C until expiration date. For more detailed information, please download the following document on our website.

Sensitivity

5 pg/mL



Tel: 1-631-624-4882 (USA)

Tel: 44-161-818-6441 (Europe)



Fax: 1-631-938-8221



Email: info@creative-diagnostics.com