



## User's Manual

# Mouse LIF ELISA Kit



DEIA2866



96T



This product is for research use only and is not intended for diagnostic use.

For illustrative purposes only. To perform the assay the instructions for use provided with the kit have to be used.

---

### Creative Diagnostics

 Address: 45-1 Ramsey Road, Shirley, NY 11967, USA

 Tel: 1-631-624-4882 (USA) 44-161-818-6441 (Europe)  Fax: 1-631-938-8221

 Email: [info@creative-diagnostics.com](mailto:info@creative-diagnostics.com)  Web: [www.creative-diagnostics.com](http://www.creative-diagnostics.com)

---

## PRODUCT INFORMATION

### Intended Use

For the quantitative determination of mouse Leukemia Inhibitory Factor (LIF) concentrations in cell culture supernates, mouse serum, and EDTA plasma.

### Storage

Store at 2-8°C. For more detailed information, please download the following document on our website.

### Precision

Intra-assay Precision:

5.1-8.5%

Inter-assay Precision:

4.3-5.7%

### Detection Range

1.34-3.13 pg/mL

### Sensitivity

1.94 pg/mL

