



User's Manual

Mouse CXCL9/MIG ELISA Kit



DEIA2882



96T



This product is for research use only and is not intended for diagnostic use.

For illustrative purposes only. To perform the assay the instructions for use provided with the kit have to be used.

Creative Diagnostics

 Address: 45-1 Ramsey Road, Shirley, NY 11967, USA

 Tel: 1-631-624-4882 (USA) 44-161-818-6441 (Europe)  Fax: 1-631-938-8221

 Email: info@creative-diagnostics.com  Web: www.creative-diagnostics.com

PRODUCT INFORMATION

Intended Use

For the quantitative determination of Mouse Monokine Induced by IFN- γ (MIG) concentrations in tissue homogenates, cell culture supernates, Mouse serum, and plasma.

Storage

Store the reagents at 2-8°C until expiration date. For more detailed information, please download the following document on our website.

Detection Range

0.98-7.8 pg/mL

Sensitivity

3.0 pg/mL

