



User's Manual

Mouse/Rat HGF ELISA Kit



DEIA2908



96T



This product is for research use only and is not intended for diagnostic use.

For illustrative purposes only. To perform the assay the instructions for use provided with the kit have to be used.

Creative Diagnostics

 Address: 45-1 Ramsey Road, Shirley, NY 11967, USA

 Tel: 1-631-624-4882 (USA) 44-161-818-6441 (Europe)  Fax: 1-631-938-8221

 Email: info@creative-diagnostics.com  Web: www.creative-diagnostics.com

PRODUCT INFORMATION

Intended Use

For the quantitative determination of mouse or rat Hepatocyte Growth Factor (HGF) concentrations in cell culture supernates, tissue homogenates, serum, and platelet-poor plasma.

Reagents And Materials Provided

1. Mouse/Rat HGF Microplate
2. Mouse/Rat HGF Conjugate
3. Mouse/Rat HGF Standard
4. Mouse/Rat HGF Control
5. Assay Diluent
6. Calibrator Diluent Concentrate
7. Wash Buffer Concentrate(25X)
8. Color Reagent A
9. Color Reagent B
10. Stop Solution

Storage

Stored at 2-8°C.

Precision

Intra-assay Precision: 3.9%-5.3%

Inter-assay Precision: 4.7%-6.8%

Detection Range

62.5-2000pg/ml

Sensitivity

Sixty-six assays were evaluated and the minimum detectable dose (MDD) of mouse/rat HGF ranged from 1.33-12.1 pg/mL. The mean MDD was 4.38 pg/mL.