



User's Manual

Mouse Fas Ligand/TNFSF6 ELISA Kit



DEIA2926



96T



This product is for research use only and is not intended for diagnostic use.

For illustrative purposes only. To perform the assay the instructions for use provided with the kit have to be used.

Creative Diagnostics

 **Address: 45-1 Ramsey Road, Shirley, NY 11967, USA**

 **Tel: 1-631-624-4882 (USA) 44-161-818-6441 (Europe)**  **Fax: 1-631-938-8221**

 **Email: info@creative-diagnostics.com**  **Web: www.creative-diagnostics.com**

PRODUCT INFORMATION

Intended Use

For the quantitative determination of mouse Fas Ligand (FasL) concentrations in cell culture supernates, mouse serum, and plasma.

Storage

Store at 4°C for frequent use, at -20°C for infrequent use. Avoid multiple freeze-thaw cycles. Return unused wells to the foil pouch containing the desiccant pack, reseal along entire edge of zip-seal. May be stored for up to 1 month at 2-8°C.

Precision

Intra-assay Precision:
4.0-9.6%

Inter-assay Precision:
4.9-8.9%

Detection Range

3.1-3.9 pg/mL

Sensitivity

3.6 pg/mL