



**User's Manual**

# Mouse Endoglin/CD105 ELISA Kit



**DEIA2944**



**96T**





This product is for research use only and is not intended for diagnostic use.

For illustrative purposes only. To perform the assay the instructions for use provided with the kit have to be used.

---

**Creative Diagnostics**

 **Address: 45-1 Ramsey Road, Shirley, NY 11967, USA**

 **Tel: 1-631-624-4882 (USA) 44-161-818-6441 (Europe)**  **Fax: 1-631-938-8221**

 **Email: [info@creative-diagnostics.com](mailto:info@creative-diagnostics.com)**  **Web: [www.creative-diagnostics.com](http://www.creative-diagnostics.com)**

---

## PRODUCT INFORMATION

### Intended Use

For the quantitative determination of mouse Endoglin concentrations in cell culture supernates, cell lysates, serum, plasma, and urine.

### Storage

Store at 2-8°C. Do not use past kit expiration date.

### Precision

Intra-assay Precision: 2.7%-5.2%

Inter-assay Precision: 4.9%-8.7%

### Detection Range

1.28-13.6 pg/mL

### Sensitivity

4.17 pg/mL

