



## User's Manual

# Human PDGF-DD ELISA Kit



DEIA2966



96T



This product is for research use only and is not intended for diagnostic use.

For illustrative purposes only. To perform the assay the instructions for use provided with the kit have to be used.

---

### Creative Diagnostics

 Address: 45-1 Ramsey Road, Shirley, NY 11967, USA

 Tel: 1-631-624-4882 (USA) 44-161-818-6441 (Europe)  Fax: 1-631-938-8221

 Email: [info@creative-diagnostics.com](mailto:info@creative-diagnostics.com)  Web: [www.creative-diagnostics.com](http://www.creative-diagnostics.com)

---

## PRODUCT INFORMATION

### Intended Use

For the quantitative determination of human Platelet-Derived Growth Factor DD (PDGF-DD) concentrations in cell culture supernates, serum, plasma, saliva, and human milk.

### Storage

Store the unopened kit at 4°C upon receipt and when it is not in use. For more detailed information, please download the following document on our website.

### Detection Range

0.699-3.51 pg/mL

### Detection Limit

1.67 pg/mL



**Tel:** 1-631-624-4882 (USA)

**Tel:** 44-161-818-6441 (Europe)



**Fax:** 1-631-938-8221



**Email:** [info@creative-diagnostics.com](mailto:info@creative-diagnostics.com)