



**User's Manual**

# Human 8-hydroxy-2'-desoxyguanosine (8-OHDG) Elisa kit

**REF**

**DEIA2984**



**96T**

**RUO**



This product is for research use only and is not intended for diagnostic use.

For illustrative purposes only. To perform the assay the instructions for use provided with the kit have to be used.

---

**Creative Diagnostics**

 **Address: 45-1 Ramsey Road, Shirley, NY 11967, USA**

 **Tel: 1-631-624-4882 (USA) 44-161-818-6441 (Europe)**  **Fax: 1-631-938-8221**

 **Email: [info@creative-diagnostics.com](mailto:info@creative-diagnostics.com)**  **Web: [www.creative-diagnostics.com](http://www.creative-diagnostics.com)**

---

## PRODUCT INFORMATION

### Intended Use

This CD 8-OHDG ELISA kit is a 1.5 hour solid-phase ELISA designed for the quantitative determination of Human 8-OHDG.

### Storage

2-8°C (6 months)

### Precision

Intra-assay Precision: 4.4%-5.6%

Inter-assay Precision: 6.6%-7.9%

### Sensitivity

The sensitivity in this assay is 1.0 ng/mL.

### Recovery

94-103%