



User's Manual

Canine SCF ELISA Kit



DEIA2987



96T



This product is for research use only and is not intended for diagnostic use.

For illustrative purposes only. To perform the assay the instructions for use provided with the kit have to be used.

Creative Diagnostics

 Address: 45-1 Ramsey Road, Shirley, NY 11967, USA

 Tel: 1-631-624-4882 (USA) 44-161-818-6441 (Europe)  Fax: 1-631-938-8221

 Email: info@creative-diagnostics.com  Web: www.creative-diagnostics.com

PRODUCT INFORMATION

Intended Use

This ELISA kit contains the basic components required for the development of sandwich ELISAs to measure natural and recombinant canine Stem Cell Factor (SCF). It is designed for the analysis of cell culture supernates. Other sample types, such as serum and plasma, need to be validated prior to use in this kit.

Storage

Store at 4°C for frequent use, at -20°C for infrequent use. Avoid multiple freeze-thaw cycles. Return unused wells to the foil pouch containing the desiccant pack, reseal along entire edge of zip-seal. May be stored for up to 1 month at 2-8°C. For more detailed information, please download the following document on our website.

