



User's Manual

Mouse IL-13 R alpha 2 ELISA Kit



DEIA2995



96T



This product is for research use only and is not intended for diagnostic use.

For illustrative purposes only. To perform the assay the instructions for use provided with the kit have to be used.

Creative Diagnostics

 **Address: 45-1 Ramsey Road, Shirley, NY 11967, USA**

 **Tel: 1-631-624-4882 (USA) 44-161-818-6441 (Europe)**  **Fax: 1-631-938-8221**

 **Email: info@creative-diagnostics.com**  **Web: www.creative-diagnostics.com**

PRODUCT INFORMATION

Intended Use

For the quantitative determination of Mouse Interleukin 13 Receptor alpha 2 (IL-13 R α 2) concentrations in cell culture supernates, tissue homogenates, serum, and plasma.

Storage

Store the kit at 2-8°C. For more detailed information, please download the following document on our website.

Detection Range

1.0-5.6 pg/mL

Sensitivity

2.3 pg/mL