



User's Manual

Mouse/Rat FABP1/L-FABP ELISA Kit

REF

DEIA2996



96T

RUO

This product is for research use only and is not intended for diagnostic use.

For illustrative purposes only. To perform the assay the instructions for use provided with the kit have to be used.

Creative Diagnostics

 **Address: 45-1 Ramsey Road, Shirley, NY 11967, USA**

 **Tel: 1-631-624-4882 (USA) 44-161-818-6441 (Europe)**  **Fax: 1-631-938-8221**

 **Email: info@creative-diagnostics.com**  **Web: www.creative-diagnostics.com**

PRODUCT INFORMATION

Intended Use

For the quantitative determination of mouse or rat Fatty Acid Binding Protein 1 (FABP1) concentrations in cell culture supernates, serum, plasma, and urine.

Storage

Store at 4°C for frequent use, at -20°C for infrequent use. Avoid multiple freeze-thaw cycles.

Return unused wells to the foil pouch containing the desiccant pack, reseal along entire edge of zip-seal. May be stored for up to 1 month at 2-8°C.

Precision

Intra-assay Precision:
2.5-3.6%

Inter-assay Precision:
3.7-4.8%

Detection Range

0.001-0.029 ng/mL

Sensitivity

0.005 ng/mL