



User's Manual

Anti-Benzo(a)pyrene Ab (Mouse) ELISA Kit

REF

DEIA3198



96T

RUO

This product is for research use only and is not intended for diagnostic use.

For illustrative purposes only. To perform the assay the instructions for use provided with the kit have to be used.

Creative Diagnostics

 **Address: 45-1 Ramsey Road, Shirley, NY 11967, USA**

 **Tel: 1-631-624-4882 (USA) 44-161-818-6441 (Europe)**  **Fax: 1-631-938-8221**

 **Email: info@creative-diagnostics.com**  **Web: www.creative-diagnostics.com**

PRODUCT INFORMATION

Intended Use

Anti-Benzo(a)pyrene Ab (Mouse) ELISA Kit is intended for qualitative or quantitative detection of anti-Benzo(a)pyrene antibodies in mouse serum.

Storage

Store the kit reagents at 2-10°C, in a dry place and protected from the light. Avoid freezing.

Precision

Intra-assay: 2.3-3.9%.

Inter-assay: 8.5-10.2%.