



User's Manual

Bisphenol A ELISA Kit



DEIA3208



96T



This product is for research use only and is not intended for diagnostic use.

For illustrative purposes only. To perform the assay the instructions for use provided with the kit have to be used.

Creative Diagnostics

 Address: 45-1 Ramsey Road, Shirley, NY 11967, USA

 Tel: 1-631-624-4882 (USA) 44-161-818-6441 (Europe)  Fax: 1-631-938-8221

 Email: info@creative-diagnostics.com  Web: www.creative-diagnostics.com

PRODUCT INFORMATION

Intended Use

Bisphenol A ELISA Kit is intended for the detection of bisphenol A (BPA) in environmental samples.

Storage

Store the kit reagents at 2-10°C, in a dry place and protected from the light. Avoid freezing. For more detailed information, please download the following document on our website.

Detection Range

0.01 µg/mL to 1 µg/mL

Detection Limit

10 ng/mL