



User's Manual

Human CA ELISA Kit



DEIA3245



96T



This product is for research use only and is not intended for diagnostic use.

For illustrative purposes only. To perform the assay the instructions for use provided with the kit have to be used.

Creative Diagnostics

 Address: 45-1 Ramsey Road, Shirley, NY 11967, USA

 Tel: 1-631-624-4882 (USA) 44-161-818-6441 (Europe)  Fax: 1-631-938-8221

 Email: info@creative-diagnostics.com  Web: www.creative-diagnostics.com

PRODUCT INFORMATION

Intended Use

This immunoassay kit allows for the in quantitative determination of Human carbonic anhydrase, CA, concentrations in serum, tissue homogenates and other biological fluids.

Storage

Store components of the kit at 2-8°C or -20°C upon arrival up to the expiration date. For more detailed information, please download the following document on our website.

Detection Range

1.56 IU/mL-100 IU /mL

Sensitivity

0.8 IU /mL

