



User's Manual

# Human EGFR (pY1086)/Pan EGFR ELISA Kit



DEIA3355



96T



This product is for research use only and is not intended for diagnostic use.

For illustrative purposes only. To perform the assay the instructions for use provided with the kit have to be used.

---

Creative Diagnostics

 Address: 45-1 Ramsey Road, Shirley, NY 11967, USA

 Tel: 1-631-624-4882 (USA) 44-161-818-6441 (Europe)  Fax: 1-631-938-8221

 Email: [info@creative-diagnostics.com](mailto:info@creative-diagnostics.com)  Web: [www.creative-diagnostics.com](http://www.creative-diagnostics.com)

---

## PRODUCT INFORMATION

### Intended Use

Human EGFR (pY1086)/Pan EGFR ELISA Kit is an enzyme-linked immunosorbent assay for the quantitative measurement of human phospho EGFR (Y1086) and Pan EGFR.

### Storage

Store the unopened kit at 4°C upon receipt and when it is not in use. For more detailed information, please download the following document on our website.

