



User's Manual

ET ELISA Kit



DEIA3405



96T



This product is for research use only and is not intended for diagnostic use.

For illustrative purposes only. To perform the assay the instructions for use provided with the kit have to be used.

Creative Diagnostics

 Address: 45-1 Ramsey Road, Shirley, NY 11967, USA

 Tel: 1-631-624-4882 (USA) 44-161-818-6441 (Europe)  Fax: 1-631-938-8221

 Email: info@creative-diagnostics.com  Web: www.creative-diagnostics.com

PRODUCT INFORMATION

Intended Use

This immunoassay kit allows for the in quantitative determination of ET concentrations in serum, plasma and other biological fluids.

Storage

The Standard, Detection Reagent A and Detection Reagent B should be stored at -20°C upon being received. Other reagents are kept according to the labels on vials. For more detailed information, please download the following document on our website.

Detection Range

0.0156-1 EU/mL

Sensitivity

0.01 EU/mL