



User's Manual

F11 (Human) ELISA Kit



DEIA3408



96T



This product is for research use only and is not intended for diagnostic use.

For illustrative purposes only. To perform the assay the instructions for use provided with the kit have to be used.

Creative Diagnostics

 Address: 45-1 Ramsey Road, Shirley, NY 11967, USA

 Tel: 1-631-624-4882 (USA) 44-161-818-6441 (Europe)  Fax: 1-631-938-8221

 Email: info@creative-diagnostics.com  Web: www.creative-diagnostics.com

PRODUCT INFORMATION

Intended Use

F11 (Human) ELISA Kit is a sandwich enzyme immunoassay for the quantitative measurement of human F11 in plasma, serum, and cell culture supernatants.

Storage

Store components of the kit at 2-8°C or -20°C upon arrival up to the expiration date. For more detailed information, please download the following document on our website.

Sensitivity

1.5 ng/mL