



User's Manual

Mouse/Rat HSV-1 IgG ELISA Kit



DEIA3555



96T



This product is for research use only and is not intended for diagnostic use.

For illustrative purposes only. To perform the assay the instructions for use provided with the kit have to be used.

Creative Diagnostics

 **Address: 45-1 Ramsey Road, Shirley, NY 11967, USA**

 **Tel: 1-631-624-4882 (USA) 44-161-818-6441 (Europe)**  **Fax: 1-631-938-8221**

 **Email: info@creative-diagnostics.com**  **Web: www.creative-diagnostics.com**

PRODUCT INFORMATION

Intended Use

HSV1 IgG (Mouse/Rat) ELISA Kit is an indirect enzyme immunoassay for the measurement of mouse/rat HSV-1 specific IgG class antibodies in serum and plasma.

Storage

Store the kit at 2-8°C. For more detailed information, please download the following document on our website.