



User's Manual

HIF1A (Human) ELISA Kit



DEIA3566



96T





This product is for research use only and is not intended for diagnostic use.

For illustrative purposes only. To perform the assay the instructions for use provided with the kit have to be used.

Creative Diagnostics

 Address: 45-1 Ramsey Road, Shirley, NY 11967, USA

 Tel: 1-631-624-4882 (USA) 44-161-818-6441 (Europe)  Fax: 1-631-938-8221

 Email: info@creative-diagnostics.com  Web: www.creative-diagnostics.com

PRODUCT INFORMATION

Intended Use

This immunoassay kit allows for the in vitro quantitative determination of human hypoxia-inducible factor 1 α , HIF-1 α concentrations in serum, plasma, urine, tissue homogenates and cell culture supernates and other biological fluids.

Storage

Store the unopened kit at 4°C upon receipt and when it is not in use. For more detailed information, please download the following document on our website.

Detection Range

0.15 ng/mL-10.0 ng/mL.

Sensitivity

0.041 ng/mL