



User's Manual

Parainfluenza 1/2/3 IgM ELISA Kit



DEIA371



96T



This product is for research use only and is not intended for diagnostic use.

For illustrative purposes only. To perform the assay the instructions for use provided with the kit have to be used.

Creative Diagnostics

 **Address:** 45-1 Ramsey Road, Shirley, NY 11967, USA

 **Tel:** 1-631-624-4882 (USA) 44-161-818-6441 (Europe)  **Fax:** 1-631-938-8221

 **Email:** info@creative-diagnostics.com  **Web:** www.creative-diagnostics.com

PRODUCT INFORMATION

Intended Use

The Parainfluenza 1/2/3 IgM antibody ELISA kit has been designed for the detection and the quantitative determination of specific IgM antibodies against Parainfluenza 1/2/3 in serum and plasma.

Storage

For more detailed information, please download the following document on our website.

Precision

Intra-assay-Precision: 5.4%

Inter-assay-Precision: 7.7%

Inter-Lot-Precision: 2.9-8.0%

Sensitivity

1.04 U/mL