



## User's Manual

# Human IL1RL1 ELISA Kit



DEIA3755



96T



This product is for research use only and is not intended for diagnostic use.

For illustrative purposes only. To perform the assay the instructions for use provided with the kit have to be used.

---

### Creative Diagnostics

 Address: 45-1 Ramsey Road, Shirley, NY 11967, USA

 Tel: 1-631-624-4882 (USA) 44-161-818-6441 (Europe)  Fax: 1-631-938-8221

 Email: [info@creative-diagnostics.com](mailto:info@creative-diagnostics.com)  Web: [www.creative-diagnostics.com](http://www.creative-diagnostics.com)

---

## PRODUCT INFORMATION

### Intended Use

The IL1RL1 (Human) ELISA (Enzyme-Linked Immunosorbent Assay) kit is an in enzyme-linked immunosorbent assay for the quantitative measurement of Human IL-1 R4 in serum, plasma, cell culture supernatants and urine.

### Storage

May be stored for up to 6 months at 2-8°C from the date of shipment. For more detailed information, please download the following document on our website.

### Sensitivity

2 pg/mL



**Tel:** 1-631-624-4882 (USA)

**Tel:** 44-161-818-6441 (Europe)



**Fax:** 1-631-938-8221



**Email:** [info@creative-diagnostics.com](mailto:info@creative-diagnostics.com)