



## User's Manual

# IL24 (Human) ELISA Kit



DEIA3771



96T



This product is for research use only and is not intended for diagnostic use.

For illustrative purposes only. To perform the assay the instructions for use provided with the kit have to be used.

---

### Creative Diagnostics

 Address: 45-1 Ramsey Road, Shirley, NY 11967, USA

 Tel: 1-631-624-4882 (USA) 44-161-818-6441 (Europe)  Fax: 1-631-938-8221

 Email: [info@creative-diagnostics.com](mailto:info@creative-diagnostics.com)  Web: [www.creative-diagnostics.com](http://www.creative-diagnostics.com)

---

## PRODUCT INFORMATION

### Intended Use

Detection and Quantification of Human Interleukin-24 (hIL-24) Concentrations in Cell Lysates, Sera and Plasma.

### Storage

Sealed, Unopened Assay Kit: 2-8°C, 1 month

Reconstituted, Opened Assay Kit: 2-8°C, 1 month For more detailed information, please download the following document on our website.

### Sensitivity

50 ng/mL