



**User's Manual**

# **Mycobacterium tuberculosis IgG ELISA Kit, Sensitive**



**DEIA384**



**96T**





This product is for research use only and is not intended for diagnostic use.

For illustrative purposes only. To perform the assay the instructions for use provided with the kit have to be used.

---

**Creative Diagnostics**

 **Address: 45-1 Ramsey Road, Shirley, NY 11967, USA**

 **Tel: 1-631-624-4882 (USA) 44-161-818-6441 (Europe)**  **Fax: 1-631-938-8221**

 **Email: [info@creative-diagnostics.com](mailto:info@creative-diagnostics.com)**  **Web: [www.creative-diagnostics.com](http://www.creative-diagnostics.com)**

---

## PRODUCT INFORMATION

### Intended Use

The Mycobacterium tuberculosis IgG antibody ELISA kit has been designed for the detection and the quantitative determination of specific IgG antibodies against Mycobacterium tuberculosis in serum and plasma.

### Storage

For more detailed information, please download the following document on our website.

### Precision

Intra-assay-Precision: 7.6%

Inter-assay-Precision: 9.4%

Inter-Lot-Precision: 3.1-9.9%

### Sensitivity

1.09 U/mL

