



User's Manual

PPBP (Human) ELISA Kit



DEIA4075



96T





This product is for research use only and is not intended for diagnostic use.

For illustrative purposes only. To perform the assay the instructions for use provided with the kit have to be used.

Creative Diagnostics

 Address: 45-1 Ramsey Road, Shirley, NY 11967, USA

 Tel: 1-631-624-4882 (USA) 44-161-818-6441 (Europe)  Fax: 1-631-938-8221

 Email: info@creative-diagnostics.com  Web: www.creative-diagnostics.com

PRODUCT INFORMATION

Intended Use

The PPBP (Human) ELISA Kit is an in vitro enzyme-linked immunosorbent assay for the quantitative measurement of human NAP-2 in serum, plasma, cell culture supernatants and urine.

Storage

All kit components of this kit are stable at 2-8°C. For more detailed information, please download the following document on our website.

Sensitivity

8.5 pg/mL



Tel: 1-631-624-4882 (USA)

Tel: 44-161-818-6441 (Europe)



Fax: 1-631-938-8221



Email: info@creative-diagnostics.com