



User's Manual

S100A7 (Human) ELISA Kit



DEIA4159



96T



This product is for research use only and is not intended for diagnostic use.

For illustrative purposes only. To perform the assay the instructions for use provided with the kit have to be used.

Creative Diagnostics

 Address: 45-1 Ramsey Road, Shirley, NY 11967, USA

 Tel: 1-631-624-4882 (USA) 44-161-818-6441 (Europe)  Fax: 1-631-938-8221

 Email: info@creative-diagnostics.com  Web: www.creative-diagnostics.com

PRODUCT INFORMATION

Intended Use

This kit is used for the quantitative measurement of human S100A7/Psoriasin in washing fluid from skin, nipple fluid, tissue culture medium and other biological samples probably including serum, plasma.

Storage

Store the unopened kit at 4°C upon receipt and when it is not in use. For more detailed information, please download the following document on our website.

Detection Range

Dilution factors need to be taken into consideration in calculating the human S100A7/Psoriasin concentration. Results exceeding human S100A7/Psoriasin level of 90 ng/mL should be repeated with diluted samples.

Detection Limit

0.12 ng/mL

