



User's Manual

SNCA (Human) ELISA Kit



DEIA4206



96T



This product is for research use only and is not intended for diagnostic use.

For illustrative purposes only. To perform the assay the instructions for use provided with the kit have to be used.

Creative Diagnostics

 Address: 45-1 Ramsey Road, Shirley, NY 11967, USA

 Tel: 1-631-624-4882 (USA) 44-161-818-6441 (Europe)  Fax: 1-631-938-8221

 Email: info@creative-diagnostics.com  Web: www.creative-diagnostics.com

PRODUCT INFORMATION

Intended Use

This immunoassay kit allows for the in vitro quantitative determination of Human α -Synuclein, SNCA, concentrations in serum, tissue homogenates and other biological fluids.

Storage

For more detailed information, please download the following document on our website.

Detection Range

0.78 ng/mL-50 ng/mL

Sensitivity

4.6 pg/mL.