



User's Manual

Varicella zoster virus IgG ELISA Kit



DEIA4425



96T



This product is for research use only and is not intended for diagnostic use.

For illustrative purposes only. To perform the assay the instructions for use provided with the kit have to be used.

Creative Diagnostics

 **Address: 45-1 Ramsey Road, Shirley, NY 11967, USA**

 **Tel: 1-631-624-4882 (USA) 44-161-818-6441 (Europe)**  **Fax: 1-631-938-8221**

 **Email: info@creative-diagnostics.com**  **Web: www.creative-diagnostics.com**

PRODUCT INFORMATION

Intended Use

Varicella zoster virus IgG ELISA assay is intended for research use of VZV associated diseases, namely varicella and herpes zoster, including their neurological complications, i.e. encephalitis, meningitis, cerebellitis, vacsulitis, myelitis or inflammatory neuropathies.

Storage

Store the kit and the kit reagents at 2-10°C, in a dry place and protected from the light. Store unused strips in the sealable pouch and keep the desiccant inside. For more detailed information, please download the following document on our website.

Sensitivity

The diagnostic sensitivity of the test is 96.6%.